

TNEED

pro



MICRONEEDLING PEN FOR MEDICAL USE ONLY

A new elegant handpiece, which works both via cable and wireless, with the possibility of using various cartridges to develop different protocols. Our PD/DS patent guarantees, as always, ease and precision of adjustment.

PRO FESSIONAL

New Motor For Even Clearer Skin Penetration



The effectiveness of a professional medical microneedling treatment lies in the creation of microchannels that enter the skin and descend into the fibroblasts.

The sharp descent of the microneedles at a controlled depth into the skin allows the production of collagen and elastin to be stimulated.

More precise is the entry and exit of the microneedles from the skin, the better is the result. **All this also allows you to reduce inflammation and pain during and after treatment.**

For this reason, 1NEED pro features a newly developed engine that is more on force and high-impact, even more professional.

PRO TOCOLS

Not Just Microneedling

Thanks to different cartridges, with 1NEED pro you can work in a combined and synergic way and carry out different treatments:

**Microneedling - Micropricking
Microdermabrasion
Specific microneedling for scars.**

PRO CEDURE

**FINE AND EXPRESSION LINES
DEEP WRINKLES · EYE CONTOUR
PIGMENTATION SPOTS AND SIGNS
OF AGE · ACNE · ROSACEA · SCARS
DILATED PORES RELAXED AND
SAGGING SKIN · STRETCH MARKS
HAIR LOSS · CELLULITIS**



PROTECTION

NOFLOW® THE SAFEST SYSTEM IN THE WORLD FOR PROTECTION AGAINST THE RISK OF CONTAMINATION



CASES BEFORE AND AFTER



**1NEED ON
HYPERPIGMENTATION**
2 NANOPEEL SESSIONS
WITH CARTRIDGE NF130AN

**1NEED ON WRINKLES
PERIOCCULAR AREA**
1 MICRONEEDLING SESSIONS
WITH CARTRIDGE NF9MN018
Working depth 0.50mm



1NEED ON STRETCH MARKS
2 MICRONEEDLING SESSIONS
WITH CARTRIDGE NF9MN018
working depth 1.75mm

1NEED ON ACNE
5 MICROPRICKING SESSIONS
WITH CARTRIDGE NF9MP018
Working depth 0.30mm



NOFLOW[®]

CARTRIDGES



World Safest Anti Back-Flow System

Cross-Contamination And Risks

The contamination of infected liquids can in fact pass from the needles which come into contact with the blood to the cartridge and from the cartridge to the handpiece.

The real and total safest solution to avoid any kind of contamination is to collect and contain any tiny infections that pass through the cartridge: NoFlow[®] Cartridges.

Thinking Out Of The Box, Or Better, Out Of The Cartridge

Existing anti back-flow solutions concentrated on creating a safety membrane after the needles inside the cartridge body. However, this solution does not guarantee absolute effectiveness because the membrane moves inside the cartridge. In this way, small infectious body fluids can still leak into the space between the membrane and the plastic body of the cartridge

towards the handpiece.

We went outside the cartridge by creating a silicone cap that collect any pathogens that may pass inside the cartridge.

In this way, at the end of the treatment, by disposing our NOFLOW[®] sterile cartridge, we will be sure that we have thrown away any possibility of cross-contamination.



International Patents:

- EU: EP3789075
- Usa: US 11,612,728 B2
- Korea 10-2023-7027166
- Usa 18/263,486
- EU 21839276.9

- **AISI304 Surgical Steel**
- **New generation of plastic polymer**
- **Strict quality production controls**
- **Batch tracking**
- **In Vitro Cytotoxicity Tests approved**
- **Hypersensitivity Test (GPMT) approved**
- **Skin Irritation test approved**
- **Validated sterilization procedure**
- **Sterile safety hard protective packaging**

A Product of

CAMPOMATS
HUMAN TECHNOLOGY

